

MULTIPLACE HYPERBARIC OXYGEN CHAMBER

OxyVista 64D

A luxurious 2-to-4 person hard-shell chamber rated to 2.0 ATA — 64" diameter × 99" length — engineered for intensive HBOT protocols.



NOISE LEVEL
< 60 dB

WARRANTY
3 Years

64"

Interior Diameter

1-4

Simultaneous Users

2.0

ATA Max Pressure

95%

O2 Delivery

State-of-the-art oxygen delivery, effortless walk-in entry, and ultimate comfort for clinics, wellness centers, and demanding home users. Redundant dual 120 LPM compressors and dual 20 LPM oxygen generators mean no external oxygen tank is ever required — just power it up, step inside, and start your session.

TRUSTED & CERTIFIED



STARTING PRICE
\$87,900

Financing available

Request a demo & quote

+1 (630) 812-7865
sales@oxygenhealthsystems.com
oxygenhealthsystems.com

KEY BENEFITS

Better outcomes for every session



Faster Recovery

Increased oxygen at 2.0 ATA promotes quicker healing for injuries, surgeries, and chronic wounds.



Improved Circulation

HBOT boosts oxygen delivery to tissues, improving circulation and cardiovascular health.



Enhanced Immunity

Higher oxygen levels strengthen the body's immune response to fight infection and recover.

ENGINEERED FEATURES

Built for safety, comfort & ease of use

Advanced AI Touch Control

Interior & exterior touchscreens show pressure, temperature, humidity & O₂.

Dual All-in-One System

Dual 120 LPM compressors & dual 20 LPM oxygen generators — no external tank.

Selectable Pressure Protocols

Incremental pressures of 1.3 / 1.6 / 1.8 / 2.0 ATA, plus optional protocols.

Easy-Slide Cabin Door

Polycarbonate sliding glass door with automatic release on de-pressurization.

Complimentary Air Conditioning

Two-part climate control with built-in dehumidifier for comfortable sessions.

Walk-In Multiplace Design

Seats 1–4 users with effortless entry and a spacious, upright interior.

OPTIONAL UPGRADES

Tailor the chamber to your practice

BIBS Oxygen Delivery

Optional

Inhale-activated 93–98% oxygen via lightweight mask; exhaled CO₂ expelled outside the chamber.

Starlight Ceiling

+\$1,200

Fiber-optic night sky with Milky Way & galaxy for a calming, spa-like session.

Premium Seating

+\$1,200/ea

Electric-reclining quilted leather chairs in gray, black, or cream finishes.

6P/OS™ Control Panel

+\$5,000

Adds six guided protocols & live safety monitoring to the standard one-touch panel.

SPECIFICATIONS

Component dimensions & power

COMPONENT	WEIGHT	POWER	DIMENSIONS
Chamber	1,720 lbs	—	99" L × 64" W × 69" H
Dual all-in-one controller	188 lbs	2,100 W	36" L × 19" W × 49" H
Air conditioning unit	60 lbs	900 W	20" L × 15" W × 23" H
Second 20 LPM O ₂ generator	170 lbs	1,760 W	36" L × 18" W × 36" H
Equipment voltage	—	110–120V*	220V on request

IDEAL FOR: Chronic wounds & ulcers · Sports injuries · Chronic pain & inflammation · Fatigue · Post-surgical recovery · Neurological conditions & TBIs

Ready to elevate your facility?

Free consultation & live demo · Custom paint & chair options · Professional installation & support

+1 (630) 812-7865

sales@oxygenhealthsystems.com
oxygenhealthsystems.com

*These statements have not been evaluated by the FDA. This product is not intended to diagnose, treat, cure, or prevent any disease.